BIOCARD™ MENOPAUSE

18.2.2002

Cat.no: 3-025-500

Rapid and reliable home test for detection of the follicle stimulating hormone (FSH) in urine.

Read the instruction for use carefully before beginning the testing.

Besides the test you need only a watch for timekeeping. Perform the test in a well-illuminated place.

INTRODUCTION


In woman's life moving from the fertile age into non-fertile age involves many changes. During this time the following symptoms can occur: period failures, irregular cycle, hot flashes, sleeplessness, headache, pain, nervousness and feeling of depression. Often reason for these symptoms is changes in hormone balance. Body has decreased the production of estrogen and progesterone and this has caused increased production of FSH (follicle stimulating hormone).

Usually the menopause will start between the ages of 46-56. The start depends on many reasons. Stress, smoking and underweight are related to early menopause. Menopause has started when the level of FSH reaches the constant and increased level > 25 mIU/ml. BIOCARD MENOPAUSE – test detects FSH hormone in urine. With this test you can find a reason for unexplained symptoms. The test result is indicative and you should consult your doctor for treatment

TEST PROCEDURE


1. Open the pouch protecting the BIOCARD Menopause Test just before use. Remove the test from the pouch.

TEST INDICATOR WINDOW


CONTROL INDICATOR WINDOW


2. Remove the cap off the sample end of the test device. Turn the test so that the perforated side of the sample end is facing you, as shown below.


3. Hold the test as shown in the drawing below and put the perforated sample end of the device under urine stream for a few seconds.

4. Withdraw the device from the urine stream and keep the sample end of the device downwards for 10-15 seconds to let extra urine drops fall down. Do not turn the sample end upwards before reading of the test results.

5. Replace the cap on the device and place the device on a flat surface with the window side up in a well-illuminated place. Read the result 10 minutes after having placed the device on a flat surface.


INTERPRETATION OF THE TEST RESULT

Read the result by comparing the intensity of the blue colour of the test and control lines (do not compare the width of the lines).


In case the test line is absent or if it is weaker (less intensely blue) than the control line, your FSH-level is not increased.

FSH LEVEL NOT INCREASED

If the test line is equally blue as the control line or if it is even more intensely coloured than the control line, your FSH-level has increased.

FSH LEVEL IS INCREASED

If you get a positive result you should repeat the testing after seven days with a new test device. If the result of the second testing is positive, the FSH level is constantly over 25 mIU/ml. If the first result is positive but a second negative, it is probable that your menopause has not yet started and the positive result was gained at the time of ovulation where there is a FSH peak.

If no lines are formed, you have likely not followed the instructions for use or the test unit has been damaged. In such a case repeat the testing with a new test unit.

WARNINGS AND LIMITATIONS

Read the instructions for use before performing a test. Failing to do so may lead to false results. Do not use expired tests or tests with a damaged pouch. Use each test device only once.

FSH concentration in the urine varies during the day. The level is at its highest in the morning. Therefore, you can obtain the most reliable result using first morning urine.

Hormonal treatment affects the FSH concentration. Hormone replacement therapy decreases the FSH level and thus may give negative results.

BIOCARD MENOPAUSE test shall be used only for *in vitro* determination of FSH level in urine according to the instruction for use. No diagnosis or decisions concerning the treatment shall be done solely on the basis of test results. If the result is contrary to your expectations, contact your gynaecologist.

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