

BIOCARD™ STREP A

A rapid, simple and reliable immunochromatographic test for the qualitative detection of group A Streptococci from a throat swab sample.

Cat. no. 3-017-400

- * Direct detection of Strep A antigens in throat swab specimens
 - * Easy and quick test procedure, result in 5 minutes
 - * Room temperature storage, shelf life 18 months
 - * Sensitivity and specificity on other side of paper
 - * **CE** – marked * 20 tests per package * Controls included
-



Helsinki, Finland
Tel. +358 - 20 155 7510
Fax. +358 - 20 155 7517
www.anibiotech.fi

BIOCARD™ STREP A

A rapid, simple and reliable immunochromatographic test for the qualitative detection of group A Streptococci from a throat swab sample.

Sensitivity:

A close-to-borderline low positive, 1.2×10^5 , gives a distinct positive result in 5 minutes.

Spesifity:

Following strains have been tested with Biocard Strep A test and found negative:

Candida albicans ATCC 28367
Cryptococcus neoformans
Escherichia coli ATCC 21432
Escherichia coli ATCC 8739
Escherichia coli ATCC 25922
Pseudomonas aeruginosa
Staphylococcus aureus ATCC 6538
Staphylococcus epidermicus
Staphylococcus saprophyticus ATCC 2866
Streptococcus agalactiae ATCC 27956
Streptococcus faecalis ATCC 19433
Streptococcus Group C

If the test can not be performed immediately the sample may be stored in a clean, dry vial for up to 72 hours. No transport media needed.



Helsinki, Finland
Tel. +358 - 20 155 7510
Fax. +358 - 20 155 7517
www.anibiotech.fi